



Pro-CT AEC

04-301



The phantom used for AEC testing on CTs. All mounting and positioning features are the same as in the Pro-CT mk II phantom for off-the-table measurement.



Technical data (can be modified to customer specifications):

- 15 PMMA ellipses 20 mm thick and with different height and width
- elliptical cross section with major and minor axes ratio equal to 3:2
- ellipse sizes range from 175,5 mm x 117 mm to 428,7 mm x 285,8 mm
- positioning stand for on and off table measurements
- carrying case

Product features:

- complies with:
 - IEC 61223-3-5
 - IEC 61223-2-6
- the Manual provides detailed guidelines for carrying out each test, results assessment and registration