

Pro-MAM AEC PMMA



- 03-201 - 180 x 240 mm version
- 03-203 - 180 x 240 mm, high precision with thickness accuracy of about 0.05 mm
- 03-204 - 240 x 300 mm version
- 03-205 - 240 x 300 mm, high precision with thickness accuracy of about 0.05 mm



Set of acrylic plates for testing Automated Exposure Control of mammography systems.

Technical data (can be modified to customer specifications):

- dimensions: 240 x 180 mm and 240 x 300 mm
- 7x 10 mm thick
- 2x 5 mm thick
- made of a transparent PMMA
- other sizes upon request

Product features:

- Complies with:
 - IEC 61223-3-2
- the Manual provides detailed guidelines for carrying out each test, results assessment and registration