

Pro-MAM Compress

03-103



Calibrated scale for assessment of the **breast compression force** of mammography systems.

Technical data (can be modified to customer specifications):

- measuring range up to 25 kg (03-103) or 34 kg (03-103-34)
- accuracy 10 g
- large, backlit, detachable (130 cm long cable) LCD display
- simple operation with three buttons
- scale size: 20 x 16 cm
- battery-operated

Product features:

- complies with:
 - IEC 61223-3-2
 - “European Guidelines for Quality Assurance in Mammography Screening”
- calibration certificate
- the Manual provides detailed guidelines for carrying out each test, results assessment and registration