


# Pro-MAM 21 Steps

03-202



**21 step aluminium wedge** for determination of **the dose reproducibility and sensitometric curve shape, speed and mid-gradient** of mammography X-ray film.

## Technical data (can be modified to customer specifications):

- dimensions: 10 x 105 x 6.3 mm
- 21 steps with a width of 5 mm
- 0.3 mm graduation per step
- made of a the highest purity aluminium
- other sizes of step wedges upon request

## Product features:

- Complies with:
  - IEC 61223-3-2
- the Manual provides detailed guidelines for carrying out each test, results assessment and registration