

Pro-MAM VOX

03-602



A breast tissue simulating test object designed to be used to compare image quality between different exposure parameters within a single system or across different breast tomosynthesis systems.

Features:

- Standard breast thickness (45 mm) 220 mm diameter semi-circle comprising breast tissue simulating material encased in PMMA.
- Microcalcifications
4 x groups of particles representing micro-calcifications randomly located in the XY plane and at approximately 10 mm intervals in the Z plane.

Please note that specific locations for calcifications cannot be given.

Pro-MAM VOX phantom (03-602, Pro-Project) allows to carry out the same tests and evaluations as VOX MAM phantom (Leeds), as per the attached statement.