

Pro-Dent 2D BASIC kit



01-001



Powered by:

DIAGNOMATIC

This kit is a versatile set of phantoms and software for carrying out constancy tests of intra-oral and OPG dental X-ray units. Thanks to the Diagnostics subscription, all tests can be quickly and effortlessly automatically analysed online.

Standard kit configuration:

- Pro-Dent set All (01-203)
- Diagnostics ECO annual subscription
- carrying case with foam inlay

The kit can be used to measure:

- dose reproducibility
- development process stability
- alignment and beam geometry
- spatial resolution
- low contrast resolution

Product features:

- Complies with:
 - OEC 61223-2-7 and IEC 61223-3-4
 - DIN 6868-5 and DIN V 6868-151
 - ÖNORM S 5240-5 and ÖNORM S 5240-11
- the manual provides detailed guidelines for carrying out each test, results assessment and registration