

# Pro-Dent All BASIC kit



01-021



Powered by:

**DIAGNOMATIC**

This kit is a versatile set of phantoms and software for carrying out constancy and acceptance tests of Intra-oral, OPG, CBCT, DVT and other 3D imaging devices. Thanks to the Diagnostics subscription all tests can be quickly and effortlessly automatically analysed online.

## Standard kit configuration:

- Pro-Dent set All (01-203)
- Pro-Dent CT MINI (01-303)
- Diagnostics BASIC annual subscription
- carrying case with dedicated foam inlay

## The kit can be used to measure:

- dose reproducibility
- development process stability
- alignment and beam geometry
- image geometry
- pixel (matrix) size
- artefacts, noise
- homogeneity
- linearity
- contrast
- low-contrast resolution (contrast sensitivity)

## Product features:

- Complies with:
  - IEC 61223-3-4, IEC 61223-3-5 and IEC 61223-2-7
  - DIN 6868-5 and DIN V 6868-151
  - ÖNORM S 5240-5 and ÖNORM S 5240-11
- the manual provides detailed guidelines for carrying out each test, results assessment and registration