



# Pro-RF BASIC kit

02-011



Powered by:



This kit is a versatile set of phantoms and software for carrying out constancy tests of radiography and fluoroscopy units. Thanks to the Pro-Control.online subscription, all tests can be quickly and effortlessly automatically analysed online and with the desktop app.

## Standard kit configuration:

- Pro-Fluo 150
- Pro-RF AI 25
- Pro-RF Rack
- Pro-RF Stand
- Pro-RF AEC Cu
- Pro-Control.online BASIC annual subscription
- carrying case with dedicated foam inlay

## The kit can be used to measure:

- collimation/beam alignment
- position and size of the effective radiation field
- dynamic range
- spatial resolution
- contrast resolution
- homogeneity
- beam quality
- image distortion
- AEC tests
- dose reproducibility

## Product features:

- complies with:
  - IEC 61223-3-1
  - DIN 6868-150
- the Manual provides detailed guidelines for carrying out each test, results assessment and registration