




# Pro-RF Stand

02-302



**A spacer plate with adjustable legs for positioning Pro-Alpha and Pro-Digi test phantoms near the X-Ray tube** on units without insertion slots (e.g. over-table / under-table tubes, C-arms etc.).

## Technical data (can be modified to customer specifications):

- table top made of PMMA
- legs made of aluminium
- dedicated space for the Pro-RF Al 25 absorption filter

## Product features:

- Complies with:
  - IEC 61223-3-1
- the Manual provides detailed guidelines for carrying out each test, results assessment and registration

