



# Pro-RF AI 21

02-307



**A 21 mm aluminum filter with insertion frame** for mounting near the X-Ray tube.

### Technical data (can be modified to customer specifications):

- filter dimensions: 120 x 120 x 21 mm
- mounting frame
- made of the purest aluminium
- optional 1 mm thick copper filter (that can be attached to the aluminium) for tests at and above 100 kV (Pro-RF Cu 1) (02-304)
- other sizes of filters upon request

### Product features:

- Complies with:
  - IEC 61223-3-1
- the Manual provides detailed guidelines for carrying out each test, results assessment and registration

