



Pro-RF MTF

02-107



The edge test device according to IEC 62220-1 **for determination of Modulation Transfer Function** that consists of a tungsten plate fixed on a lead plate.

Technical data (can be modified to customer specifications):

- overall dimensions: 220 x 110 x 13 mm
- tungsten plate dimensions: 100 x 80 x 1 mm
- lead plate dimensions: 200 x 100 x 3 mm
- tungsten edge scragginess less than 5 μm
- PMMA cover
- three leveling screws

Product features:

- Complies with:
 - IEC 62220-1
 - DIN EN 62220-1
- the Manual provides detailed guidelines for carrying out each test, results assessment and registration