




Pro-RF Ruler

- 02-601 - PMMA ruler with the radiopaque scale - 50 cm / 2 mm scale
- 02-602 - PMMA ruler with the radiopaque scale - 50 cm / 1 mm scale
- 02-603 - aluminium ruler with the radiopaque scale - 50 cm / 2 mm scale
- 02-604 - aluminium ruler with the radiopaque scale - 50 cm / 1 mm scale
- 02-605 - two PMMA rulers with the radiopaque scale joined in the middle (cross)- 40 cm each / 1 mm scale



The radiopaque ruler in the image provides a permanent reference of image size . It can be used to check the accuracy and position of the light field for quality control and adjustment as well as be included during x-ray procedures, leaving a clear measurement right on the x-ray image itself.

Technical data (can be modified to customer specifications):

We manufacture many sizes and types of radiopaque rulers for many applications, with different scales and accuracy. Please contact us to find the best solution for you.

Product features:

- Complies with:
 - IEC 61223-3-1
- the Manual provides detailed guidelines for carrying out each test, results assessment and registration