



# Pro-RF Tomo

02-116



**A set of simple tools dedicated for measurements on tomography X-Ray equipment** . It can be used for the layer height adjustment and tomography movement tests.

## Technical data (can be modified to customer specifications):

- 160 mm ruler made of PMMA with engraved detailed X-Ray opaque scale
- ruler can be folded to a very compact size for easy transportation
- stainless steel step with holes for displaying tomographic movement

## Product features:

- Complies with:
  - IEC 61223-3-1
- the Manual provides detailed guidelines for carrying out each test, results assessment and registration