



# Pro-RF Res 150

02-120



The phantom for acceptance and constancy testing of **Cone-Beam CT / Volume Tomography (DVT) in fluoroscopy** according to **DIN 6868-150**

## Technical data (can be modified to customer specifications):

- dimensions: 120 x 120 mm
- 10mm thick test module containing sets of cylindrical holes, parallel to Z-Axis, of 0.5, 0.6, 0.7, 0.8, 0.9, 1.0 and 1.3 mm diameters
- two additional 25 mm PMMA plates
- carrying case

## Product features:

- complies with:
  - DIN 6868 – 150
- the Manual provides detailed guidelines for carrying out each test, results assessment and registration