



# Pro-RF AAPM 31

02-205



**A set of acrylic and aluminium filters** for testing diagnostic X-rays according to **AAPM recommendations from Report No. 31** "Standardized Methods for Measuring Diagnostic X-Ray Exposure". Patient-equivalent acrylic and aluminium phantoms provide the necessary attenuation between the source and AEC or ABC detectors.

## Technical data (can be modified to customer specifications):

- 5x acrylic plate 250 x 250 x 25,4 mm
- 1x acrylic plate 250 x 250 x 50,8 mm
- 1x 250 x 250 x 1 mm aluminium filter
- 1x 250 x 250 x 2 mm aluminium filter
- 1x 70 x 250 x 4,5 mm aluminium filter
- 4x positioning rod
- 4x 50,8 mm spacer tube
- optional heavy-duty transport case (02-206)

## Product features:

- complies with:
  - report No. 31 "Standardized Methods for Measuring Diagnostic X-Ray Exposure"
- the Manual provides detailed guidelines for carrying out each test, results assessment and registration