



Pro-RF AEC PMMA

02-201 - 240 x 240 mm version
02-221 - 300 x 300 mm version



A set of acrylic plates for testing Automated Exposure Control of radiography equipment.

Technical data (can be modified to customer specifications):

- dimensions: 240 x 240 mm
- 9x 20 mm thick
- 1x 10 mm thick
- 2x 5 mm thick
- made of transparent PMMA
- optional high purity 240 x 240 x 25 mm aluminium plate (02-202)
- other sizes of filters upon request

Product features:

- Complies with:
 - IEC 61223-3-1
- the Manual provides detailed guidelines for carrying out each test, results assessment and registration