



Pro-RF AEC Cu

02-203



A set of high purity copper plates for testing Automated Exposure Control of radiography equipment. **It can be mounted near the X-Ray tube** and used instead of PMMA plates.

Technical data (can be modified to customer specifications):

- dimensions: 169 x 176 mm
- 1x 1.0 mm thick
- 1x 0.5 mm thick
- made of high purity copper
- other sizes of filters upon request

Product features:

- Complies with:
 - IEC 61223-3-1
- the Manual provides detailed guidelines for carrying out each test, results assessment and registration