



Pro-RF FocalSpot

02-501



A simple test device for the focal spot size evaluation . It consists of a precision bar pattern placed at the end of the cylinder that provides accurate and reproducible target-to-image receptor spacing.

Technical data (can be modified to customer specifications):

- cylinder dimensions: 152,4 x Ø76 mm
- bar pattern mounted on the top of the cylinder
- 8 resolution bar patterns ranging from 0,84 to 5,66 LP/mm

Product features:

- Complies with:
 - IEC 61223-3-1
- the Manual provides detailed guidelines for carrying out each test, results assessment and registration