



Pro-RF 21 Steps

02-305



A 21 step aluminium wedge for determination of the **dose reproducibility and sensitometric curve shape, speed and mid-gradient** of X-Ray screen-film radiography systems.

Technical data (can be modified to customer specifications):

- dimensions: 231 x 110 x 31.5 mm
- 21 steps with a width of 11 mm
- 1.5 mm graduation per step
- made of a highest purity aluminium
- copper numbers marking each step
- other sizes of step wedges upon request

Product features:

- Complies with:
 - IEC 61223-3-1
- the Manual provides detailed guidelines for carrying out each test, results assessment and registration