



# Pro-Slit

05-101



The Pro-Slit phantom is the so called **slit camera for accurate measurement of the focal spot size according to IEC 60336:2005**. Its main advantages are **repeatability, accuracy and possibility to measure the size of any focal spot**. Because of the small attenuation of the outer cover this tool can also be used with dental digital detectors - detector receives enough radiation to be triggered.

## Technical data (can be modified to customer specifications):

- dimensions: 40 x 40 x 11.5 mm
- 0.01 mm wide and 10 mm long slit with 8° spread
- diaphragm made of tungsten
- embedded in PMMA with additional stainless steel shielding
- optional Pro-Stand in a heavy duty carrying case with two layer foam inly

## Product features:

- complies with:
  - IEC 60336:2005
  - IEC 61223-3-1, 2, 4
- the Manual provides detailed guidelines for carrying out each test, results assessment and registration