



# Pro-Stand

05-102 - R/F applications

05-103 - All applications including mammography



**Adjustable stand**, designed to make focal spot measuring procedures **easy to perform**, as well as **ensuring accurate results**. It can also be used for the HVL measurements.

## Technical data (can be modified to customer specifications):

- adjustable height from 350 mm to over 600 mm (wide range of magnification)
- adjustable horizontal position and vertical angle (optional - for testing mammography systems) of the table top
- place in the base for the Pro-Dent positioning stand for accurate positioning of dental RVG detectors
- optional add-ons for stand alignment verification – especially useful for measurements on mammography units
- optional wider arm for Pro-HVL filters
- optional heavy duty carrying case with two layer foam inly

## Product features:

- complies with:
  - IEC 60336:2005
  - IEC 61223-3-1, 2, 4
- the Manual provides detailed guidelines for carrying out each test, results assessment and registration