




Pro-HVL

05-201 - standard set
05-202 - mammography set



Different **sets of high (99.9%) and standard (99.5%) purity aluminium plates** for testing **Half Value Layer**

Technical data (can be modified to customer specifications):

- dimensions: 89 x 89 mm
- mammography set: 6 x 0.1 mm, aluminium purity 99.9%
- standard set: 5 x 1.0 mm, 2 x 0.5 mm, 4 x 0.1 mm, aluminium purity 99.5%
- other sizes upon request

Product features:

- complies with:
 - IEC 61223-3-1, 2, 4
- the Manual provides detailed guidelines for carrying out each test, results assessment and registration