



Pro-HVL Cu

05-203



Different sets of **high purity copper plates** for testing **Half Value Layer**

Technical data (can be modified to customer specifications):

- dimensions: 100 x 100 mm
- standard set: 1x 2.0 mm, 2x 1.0 mm, 1x 0.5 mm, 1x 0.25 mm, 4x 0.1 mm
- other sizes upon request

Product features:

- complies with:
 - IEC 61223-3-1, 2, 4
- the Manual provides detailed guidelines for carrying out each test, results assessment and registration